

Blunt Fill Needle

Blunt Fill Needle with 5 Micron Filter

The Blunt Fill Needle and Blunt Fill Needle with integral 5 micron filter are used for draw-up of medication, re-constitution of medication and accessing IV bags. The integral 5 micron filter membrane is designed to prevent larger particles from being drawn into the syringe when preparing medications from ampoules.

Features and Benefits

45° bevel for ease in draw-up of medication from vial

18G for ease in draw-up of viscous medications

Conventional standard hub to fit standard luer lock or luer slip syringes

The BN1815F has an integral 5 micron filter for effective filtration of contaminant particles

Color coded for identification. Helps lower costs versus sharps engineered devices

Packaging

BN1815 Blunt Fill Needle	SHIPPER	1000 units per case (545mm x 188mm x 116mm)
	SHELF PACK	100 each (105mm x 84mm x 86mm)
	EXPIRATION	5 Years
BN1815F Blunt Fill Needle with 5 Micron Filter	SHIPPER	1000 units per case (545mm x 188mm x 125mm)
	SHELF PACK	100 each (105mm x 84mm x 93mm)
	EXPIRATION	5 years



PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Distributed by:

Smiths Medical International
1500 Eureka Park, Lower Pemberton
Ashford Kent, TN25 4BF
Tel: +44 (0)845 850 0445
www.smiths-medical.com

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The Smiths Medical design mark is a trademark of Smiths Medical. ©2016 Smiths Medical. All rights reserved. SS195236GB-072016

MKEECA-0098

CE Rx
ONLY

smiths medical