

smiths medical
portex®

Portex® Arterial Blood Sampling Kits

Safe Application. Simple Activation. Consistent Results.



PORTEX® ARTERIAL BLOOD SAMPLING KITS



SIMPLE ACTIVATION
OF SAFETY
DEVICE REDUCES
POSSIBILITY OF
NEEDLESTICK
INJURIES

Pro-Vent™ Arterial Blood Sampling Kit

Large venting filter media for quick filling to minimize patient discomfort

NEEDLE PRO

- Engages the needle with a simple one-handed technique, freeing the other hand to attend to the site and patient.
- A single safety activation mechanism to help prevent needlestick injuries.

ADVANCED FORMULA HEPARIN

- Available with either 1 mL or 3 mL syringes and with either advanced formula dry lithium heparin or low heparin formula.
- Allows for, consistent and accurate analyzer readings, while reducing the potential for clots.

PRO-VENT® BLOOD SAMPLING KIT ADVANTAGES

- Fully depressed syringe plunger position provides ease of aspiration in code situations.
- Positive stop on 3 mL syringe eliminates plunger disengagement.
- All Pro-Vent® kits with pre-attached needles are packaged in sterile pouches for use in a sterile field.



More accurate analyzer readings

Safer sample transport

Minimizes patient discomfort



CONSISTENT
RESULTS WITH
UNIQUE HEPARIN
FORMULATION FOR
OVER 20 YEARS AS
INDUSTRY STANDARD



Portex® Line Draw Plus Syringes Filter-Pro® air bubble removal device

Needleless blood sampling directly from an arterial line

DRY HEPARIN FORMULATION

- Allows for consistent and accurate readings.
- Available with either 1 mL Luer slip, 3 mL Luer lock, or 3 mL Luer slip syringes.

LOW HEPARIN FORMULATION

- Allows the full range of blood gas, electrolyte, and metabolite test results from very small blood samples without heparin interference.
- Available with either 1 mL Luer slip, 3 mL Luer lock, or 3 mL Luer slip syringe.

FILTER-PRO® DEVICE

- Safely expels air bubbles from the sample without exposure to blood or bloodborne pathogens.
- Caps blood-filled Luer lock and Luer slip syringes for added safety during transport.

Blood sampling directly from an arterial line

Easy connection to the stopcock

Ideal for use in the Labor and Delivery and NICU



ORDERING INFORMATION



Pro-Vent® Plus Arterial Blood Sampling Kits						
Product Code	Description	Needle	Volume	Heparin	Units per case	ARTG
4558PE	Pro-Vent® Plus 3mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	23G x 1.00"	3 mL	23.5 IU/mL	200	140416
4558PG	Pro-Vent® Plus 3mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	23G x 1.00"	3 mL	23.5 IU/mL	200	140416
4590PG	Pro-Vent® Plus 3mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	22G x 1.00"	3 mL	23.5 IU/mL	200	140416
4616E	Pro-Vent® Plus 1mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	-	1 mL	23.5 IU/mL	200	140416
4618PE	Pro-Vent® Plus 1mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	25G x 5/8"	1 mL	23.5 IU/mL	200	140416
4619LH	Pro-Vent® 1mL Arterial Blood Sampling Kit, Drylith Safety Device	-	1 mL	23.5 IU/mL	200	140416
4619PE	Pro-Vent® Plus 1mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	25G x 5/8"	1 mL	23.5 IU/mL	200	140416
4628PE	Pro-Vent® Plus 1mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	23G x 1.00"	1 mL	23.5 IU/mL	200	140416
4629PE	Pro-Vent® Plus 1mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	23G x 1.00"	1 mL	23.5 IU/mL	200	140416
4645PE	Pro-Vent® Plus 3mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	22G x 1.00"	3 mL	23.5 IU/mL	200	140416
4646E	Pro-Vent® Plus 3mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	-	3 mL	23.5 IU/mL	200	140416
4648PE	Pro-Vent® Plus 3mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	22G x 1.00"	3 mL	23.5 IU/mL	200	140416
4658PE	Pro-Vent® Plus 3mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	23G x 1.00"	3 mL	23.5 IU/mL	200	140416
4698LH	Pro-Vent® 3mL Arterial Blood Sampling Kit, Drylith Safety Device	-	3 mL	23.5 IU/mL	200	140416
4698PE	Pro-Vent® Plus 3mL Arterial Blood Sampling Kit, NeedlePro® Safety Device	23G x 1.00"	3 mL	23.5 IU/mL	200	140416



Line Draw Plus Arterial Blood Sampling Kits						
Product Code	Description	Needle	Volume	Heparin	Units per case	ARTG
4042LH	Line Draw Plus Arterial Blood Sampling Kit, 3mL Luer Slip	-	3 mL	23.5 IU/mL	200	140416

For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

MMSPCA-0589

Smiths Medical Australasia Pty Ltd
Suite 2.03, 97 Waterloo Road
Macquarie Park NSW 2113
www.smiths-medical.com

Distributed by:
USL Medical
494 Rodebank Road
Avondale, 1026 New Zealand
Freecall: 0800 658 814

**Rx
ONLY**
smiths medical

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Pro-Vent, Needle-Pro, Filter-Pro, Point-Lok, and the Portex and Smiths medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. Product referenced is CE-marked. ©2018 Smiths Medical. All rights reserved. SS195281GBAUS-122018