

smiths medical
jelco®




Jelco® I.V. Catheter

The trusted name in Peripheral
IV catheters globally



Jelco® I.V. Catheter



Sharp, beveled needle tip delivers consistent, easy insertions

- Thinwall catheter design provides high flow rates allowing for the use of a smaller gauge device
- FEP polymer catheter and V-point needle design allows reduction of insertion force which minimizes patient discomfort
- Consistent trim accuracy helps to improve vein access

Trusted by healthcare professionals for over 30 years

- Jelco® IV Catheters are used in pediatric/neonates for the high flow rates
- ORs use Jelco® IV Catheters for their quick flash, high flow rates, and ease of use for arterial lines and for I.V. access
- A high-performance product from the brand you trust, Jelco®

| Jelco IV Catheter | Radiopaque FEP polymer, straight hub | Clear, non-radiopaque FEP polymer, straight hub |
|---------------------------|--------------------------------------|---|
| Catheter Gauge and Length | Product Code | Product Code |
| 14G x 50mm | 4038 | 4018 |
| 16G x 50mm | 4032 | 4012 |
| 18G x 45mm | 4034 | 4014 |
| 18G x 32mm | 4035 | - |
| 20G x 45mm | 4039 | - |
| 20G x 32mm | 4036 | 4016 |
| 22G x 25mm | 4030 | 4010 |
| 24G x 19mm | 4033 | 4013 |

For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical ASD, Inc.
 6000 Nathan Lane North
 Minneapolis, MN, 55442, USA
 Tel: 1-614-210-7300
 Toll-Free USA 1-800-258-5361
www.smiths-medical.com

MS0UCA-1111

smiths medical

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Jelco and the Smiths Medical and design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2018 Smiths Medical. All rights reserved. SS195312GBSEA-032018