

smiths medical

bringing technology to life



DeltaVen® Closed System Catheter

Feel the difference, see the difference™



FEEL THE DIFFERENCE, SEE THE DIFFERENCE™


FEEL THE INSERTION OF A CLOSED SYSTEM CATHETER WITH EASY THREADING.

SEE THE FULL RANGE OF CONFIGURATIONS AND GAUGE SIZES TO MEET YOUR NEEDS.



The Right Sized Catheter For Your Patient

DeltaVen® helps meet INS recommendations to “select the smallest gauge² peripheral catheter that will accommodate the prescribed therapy and patient need” with sizes ranging from 16-26G and lengths from 3/4” to 1 3/4”.



THREE CONFIGURATIONS TO MEET YOUR HOSPITAL'S CURRENT GUIDELINES



Single Port with End Caps



Dual Port with End Caps














Dual Port with Dual Needleless Connector

REFERENCES:












1. Gonzalez Lopez JL, et al., Indwell times, complications and costs of open vs closed safety peripheral intravenous catheters: a randomized study, Journal of Hospital Infection (2013), <http://dx.doi.org/10.1016/j.jhin.2013.10.008>
2. INS Standards. Vol. 39, #1S, P.S51, Practice Criteria, 1.B

ORDERING INFORMATION












Single Port with End Caps - 100 units per case

Product Code	Gauge	Catheter Length (in)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3804272	26 	3/4	0.2	15	N/A
3814272	24 	3/4	0.23	24	330
3824272	22 	3/4	0.26	37	330
3825272	22 	1	0.27	31	330
3835272	20 	1	0.27	47	330
3836272	20 	1 1/4	0.29	42	330
3837272	20 	1 3/4	0.3	38	330
3845272	18 	1	0.3	65	330
3846272	18 	1 1/4	0.3	63	330
3847272	18 	1 3/4	0.32	58	330
3856272	16 	1 1/4	0.42	125	330

Dual Port with End Caps - 100 units per case

Product Code	Gauge	Catheter Length (in)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3800072	26 	3/4	0.25	15	N/A
3810072	24 	3/4	0.3	24	330
3820072	22 	3/4	0.3	37	330
3821072	22 	1	0.31	31	330
3831072	20 	1	0.32	47	330
3832072	20 	1 1/4	0.32	42	330
3833072	20 	1 3/4	0.33	38	330
3841072	18 	1	0.35	65	330
3842072	18 	1 1/4	0.35	63	330
3843072	18 	1 3/4	0.37	58	330
3852072	16 	1 1/4	0.47	125	330

Dual Port with Dual Needleless Connector - 80 units per case

Product Code	Gauge	Catheter Length (in)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3800572	26 	3/4	0.29	15	N/A
3810572	24 	3/4	0.34	24	330
3820572	22 	3/4	0.34	37	330
3821572	22 	1	0.35	31	330
3831572	20 	1	0.37	47	330
3832572	20 	1 1/4	0.37	42	330
3833572	20 	1 3/4	0.38	38	330
3841572	18 	1	0.39	65	330
3842572	18 	1 1/4	0.39	63	330
3843572	18 	1 3/4	0.41	58	330
3852572	16 	1 1/4	0.51	125	330

* All accessories (stopcock and valve) must be disconnected in order to achieve a pressure of 330 psi

For more information visit our website at www.smiths-medical.com/deltaven

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Distributed By:
Smiths Medical ASD, Inc.
 6000 Nathan Lane North
 Minneapolis, MN 55442, USA
 Tel: 1-614-210-7300
 Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

Manufactured by:
Delta Med SpA
 Via G. Rossa, 20-46019
 Viadana (Mn) - Italy
 0375 785915

MMSPCA-0014

**Rx
 ONLY**

smiths medical

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2018 Smiths Medical. All rights reserved. SS195451EN-032018