

**smiths medical**

bringing technology to life



## DeltaVen<sup>®</sup> Closed System Catheter

Feel the difference, see the difference<sup>™</sup>



# FEEL THE DIFFERENCE, SEE THE DIFFERENCE™


FEEL THE INSERTION OF A CLOSED SYSTEM CATHETER WITH EASY THREADING.

SEE THE FULL RANGE OF CONFIGURATIONS AND GAUGE SIZES TO MEET YOUR NEEDS.



## The Right Sized Catheter For Your Patient

DeltaVen® helps meet INS recommendations to “select the smallest gauge<sup>2</sup> peripheral catheter that will accommodate the prescribed therapy and patient need” with sizes ranging from 16-26G and lengths from 3/4” to 1 3/4”



### THREE CONFIGURATIONS TO MEET YOUR HOSPITAL'S CURRENT GUIDELINES



Single Port with End Caps



Dual Port with End Caps














Dual Port with Dual Needleless Connector

#### REFERENCES:












1. Gonzalez Lopez JL, et al., Indwell times, complications and costs of open vs closed safety peripheral intravenous catheters: a randomized study, Journal of Hospital Infection (2013), <http://dx.doi.org/10.1016/j.jhin.2013.10.008>
2. INS Standards. Vol. 39, #1S, P.S51, Practice Criteria, 1.B

## ORDERING INFORMATION












### Single Port with End Caps - 25 units per box, 100 units per case

Product Code	Gauge	Catheter Length (in)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3804222	26 	3/4	0.2	15	N/A
3814222	24 	3/4	0.23	24	330
3824222	22 	3/4	0.26	37	330
3825222	22 	1	0.27	31	330
3835222	20 	1	0.27	47	330
3836222	20 	1 1/4	0.29	42	330
3837222	20 	1 3/4	0.3	38	330
3845222	18 	1	0.3	65	330
3846222	18 	1 1/4	0.3	63	330
3847222	18 	1 3/4	0.32	58	330
3856222	16 	1 1/4	0.42	125	330

### Dual Port with End Caps - 25 units per box, 100 units per case

Product Code	Gauge	Catheter Length (in)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3800022	26 	3/4	0.25	15	N/A
3810022	24 	3/4	0.3	24	330
3820022	22 	3/4	0.3	37	330
3821022	22 	1	0.31	31	330
3831022	20 	1	0.32	47	330
3832022	20 	1 1/4	0.32	42	330
3833022	20 	1 3/4	0.33	38	330
3841022	18 	1	0.35	65	330
3842022	18 	1 1/4	0.35	63	330
3843022	18 	1 3/4	0.37	58	330
3852022	16 	1 1/4	0.47	125	330

### Dual Port with Dual Needleless Connector - 20 units per box, 80 units per case

Product Code	Gauge	Catheter Length (in)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3800522	26 	3/4	0.29	15	N/A
3810522	24 	3/4	0.34	24	330
3820522	22 	3/4	0.34	37	330
3821522	22 	1	0.35	31	330
3831522	20 	1	0.37	47	330
3832522	20 	1 1/4	0.37	42	330
3833522	20 	1 3/4	0.38	38	330
3841522	18 	1	0.39	65	330
3842522	18 	1 1/4	0.39	63	330
3843522	18 	1 3/4	0.41	58	330
3852522	16 	1 1/4	0.51	125	330

\* All accessories (stopcock and valve) must be disconnected in order to achieve a pressure of 330 psi

For more information visit our website at [www.smiths-medical.com/deltaven](http://www.smiths-medical.com/deltaven)

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN OTHER COUNTRIES

**Distributed By:**  
**Smiths Medical ASD, Inc.**  
 6000 Nathan Lane North  
 Minneapolis, MN 55442, USA  
 Tel: 1-614-210-7300  
 Toll-Free USA: 1-800-258-5361  
[www.smiths-medical.com](http://www.smiths-medical.com)

**Manufactured by:**  
**Delta Med SpA**  
 Via G. Rossa, 20-46019  
 Viadana (Mn) - Italy  
 0375 785915

MMSPCA-0237

**Rx  
 ONLY**

smiths medical

Find your local contact information at: [www.smiths-medical.com/customer-support/contact-us](http://www.smiths-medical.com/customer-support/contact-us)

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2018 Smiths Medical. All rights reserved. SS195451ENCA-042018