

smiths medical

bringing technology to life



DeltaVen® Closed System Catheter

Feel the difference, see the difference™



FEEL THE DIFFERENCE, SEE THE DIFFERENCE™

FEEL THE INSERTION OF A CLOSED SYSTEM CATHETER WITH EASY THREADING.


SEE THE FULL RANGE OF CONFIGURATIONS AND GAUGE SIZES TO MEET YOUR NEEDS.



Featuring a 26G Closed System Catheter

IDEAL FOR NEONATES, PEDIATRICS AND ONCOLOGY

Be confident with a range of catheters designed with the smallest patients in mind. The DeltaVen® closed system catheter is the only catheter designed to meet the unique needs of the smallest veins when critical care is needed.



FOUR CONFIGURATIONS TO BEST MEET YOUR HOSPITAL'S CURRENT GUIDELINES



Single Port



Single Port with Stopcock and Needleless Connector



Dual Port with Dual End Caps



Dual Port with Dual Needleless Connector

REFERENCES:

1. Gonzalez Lopez JL, et al., Indwell times, complications and costs of open vs closed safety peripheral intravenous catheters: a randomized study, *Journal of Hospital Infection* (2013), <http://dx.doi.org/10.1016/j.jhin.2013.10.008>.

ORDERING INFORMATION

Single Port - 25 units per case, 100 units per carton

Product Code	Gauge	Catheter Length (mm)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3804222	26	19	0.2	15	N/A
3814222	24	19	0.23	24	330
3824222	22	19	0.26	37	330
3825222	22	25	0.27	31	330
3835222	20	25	0.27	47	330
3836222	20	32	0.29	42	330
3837222	20	45	0.3	38	330
3845222	18	25	0.3	65	330
3846222	18	32	0.3	63	330
3847222	18	45	0.32	58	330
3856222	16	32	0.42	125	330

Single Port with Stopcock and Needleless Connector - 20 units per case, 80 units per carton

Product Code	Gauge	Catheter Length (mm)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3804422	26	19	0,32	15	N/A
3814422	24	19	0,37	24	330
3824422	22	19	0,38	37	330
3825422	22	25	0,38	31	330
3835422	20	25	0,39	47	330
3836422	20	32	0,39	42	330
3837422	20	45	0,41	38	330
3845422	18	25	0,42	65	330
3846422	18	32	0,42	63	330
3847422	18	45	0,44	58	330
3856422	16	32	0,54	125	330

Dual Port with End Caps - 25 units per case, 100 units per carton

Product Code	Gauge	Catheter Length (mm)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3800022	26	19	0.25	15	N/A
3810022	24	19	0.3	24	330
3820022	22	19	0.3	37	330
3821022	22	25	0.31	31	330
3831022	20	25	0.32	47	330
3832022	20	32	0.32	42	330
3833022	20	45	0.33	38	330
3841022	18	25	0.35	65	330
3842022	18	32	0.35	63	330
3843022	18	45	0.37	58	330
3852022	16	32	0.47	125	330

Dual Port with Dual Needleless Connectors - 20 units per case, 80 units per carton

Product Code	Gauge	Catheter Length (mm)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3800522	26	19	0.29	15	N/A
3810522	24	19	0.34	24	330
3820522	22	19	0.34	37	330
3821522	22	25	0.35	31	330
3831522	20	25	0.37	47	330
3832522	20	32	0.37	42	330
3833522	20	45	0.38	38	330
3841522	18	25	0.39	65	330
3842522	18	32	0.39	63	330
3843522	18	45	0.41	58	330
3852522	16	32	0.51	125	330

* All accessories (stopcock and valve) must be disconnected in order to achieve a pressure of 330 psi

For more information visit our website at www.smithsmedical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Distributed By:
Smiths Medical International
 1500 Eureka Park, Lower Pemberton
 Ashford Kent, TN25 4BF
 Tel: +44 (0)845 850 0445
www.smiths-medical.com

Manufactured by:
Delta Med SpA
 Via G. Rossa, 20-46019
 Viadana (Mn) - Italy
 0375 785915

MPAUC-1675

**Rx
 ONLY**

smiths medical

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2017 Smiths Medical. All rights reserved. SS195451GB-112017