

smiths medical

bringing technology to life



DeltaVen® Closed System Catheter

Feel the difference, see the difference™



FEEL THE DIFFERENCE, SEE THE DIFFERENCE™

FEEL THE INSERTION OF A CLOSED SYSTEM CATHETER WITH EASY THREADING.


SEE THE FULL RANGE OF CONFIGURATIONS AND GAUGE SIZES TO MEET YOUR NEEDS.



Featuring a 26G Closed System Catheter

IDEAL FOR NEONATES, PEDIATRICS AND ONCOLOGY

Be confident with a range of catheters designed with the smallest patients in mind. The DeltaVen® closed system catheter is the only catheter designed to meet the unique needs of the smallest veins when critical care is needed.



THREE CONFIGURATIONS TO BEST MEET YOUR HOSPITAL'S CURRENT GUIDELINES



Single Port



Dual Port with Dual End Caps














Dual Port with Single Needleless Connector

REFERENCES:












1. Gonzalez Lopez JL, et al., Indwell times, complications and costs of open vs closed safety peripheral intravenous catheters: a randomized study, Journal of Hospital Infection (2013), <http://dx.doi.org/10.1016/j.jhin.2013.10.008>.

ORDERING INFORMATION












Single Port - 25 units per case, 100 units per carton

Product Code	Gauge	Catheter Length (mm)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3804222	26 	19	0.2	15	N/A
3814222	24 	19	0.23	24	330
3824222	22 	19	0.26	37	330
3825222	22 	25	0.27	31	330
3835222	20 	25	0.27	47	330
3836222	20 	32	0.29	42	330
3837222	20 	45	0.3	38	330
3845222	18 	25	0.3	65	330
3846222	18 	32	0.3	63	330
3847222	18 	45	0.32	58	330
3856222	16 	32	0.42	125	330

Dual Port with Dual End Caps - 25 units per case, 100 units per carton

Product Code	Gauge	Catheter Length (mm)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3800022	26 	19	0.25	15	N/A
3810022	24 	19	0.3	24	330
3820022	22 	19	0.3	37	330
3821022	22 	25	0.31	31	330
3831022	20 	25	0.32	47	330
3832022	20 	32	0.32	42	330
3833022	20 	45	0.33	38	330
3841022	18 	25	0.35	65	330
3842022	18 	32	0.35	63	330
3843022	18 	45	0.37	58	330
3852022	16 	32	0.47	125	330

Dual Port with Single Needleless Connector - 20 units per case, 80 units per carton

Product Code	Gauge	Catheter Length (mm)	Priming Volume (mL)	Flow (mL/min)	PSI Rating
3800122	26 	19	0.29	15	N/A
3810122	24 	19	0.34	24	330
3820122	22 	19	0.34	37	330
3821122	22 	25	0.35	31	330
3831122	20 	25	0.37	47	330
3832122	20 	32	0.37	42	330
3833122	20 	45	0.38	38	330
3841122	18 	25	0.39	65	330
3842122	18 	32	0.39	63	330
3843122	18 	45	0.41	58	330
3852122	16 	32	0.51	125	330

* All accessories (stopcock and valve) must be disconnected in order to achieve a pressure of 330 psi

For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Distributed By:

Smiths Medical Australasia Pty. Ltd.

Suite 201, 10 Norbrik Drive
Bella Vista, NSW 2153
Tel: AU: 1-800-654-949
NZ: 0-800-444-200

www.smiths-medical.com

USL Medical

494 Rodebank Road
Avondale, 1026 New Zealand
Freecall: 0800 658 814

Manufactured by:

Delta Med SpA

Via G. Rossa, 20-46019
Viadana (Mn) - Italy
0375 785915

MMSPCA-0106

Rx CE
ONLY 0123

smiths medical

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2018 Smiths Medical. All rights reserved. SS195451GBAUS-022018