

smiths medical
bringing technology to life



DeltaVen® Closed System Catheter

Feel the difference,
see the difference™



FEEL THE DIFFERENCE, SEE THE DIFFERENCE™

A range of catheters designed with the smallest patient in mind. The DeltaVen® closed system catheter was created to meet the unique needs of these patients when critical care is needed.

FEEL THE DIFFERENCE FOR THE SMALLEST PATIENTS

- Taking care of the most delicate patients requires precision and the right tools. For neonates and patients with fragile veins, catheter size shouldn't be an issue. That's why we created the DeltaVen® 26 gauge safety catheter.
- With neonates, every drop of fluid counts. The 26G catheter has priming volume as low as 0.2 mL in the single port configuration.¹

SEE THE DIFFERENCE IN PATIENT CARE

- A 26G safety catheter eliminates the choice between safety and catheter size.
- DeltaVen® helps you meet INS recommendations to "select the smallest gauge peripheral catheter that will accommodate the prescribed therapy and patient need ²."
- Bonded extension set reduces the risk of contamination from repeated disconnections and reconnections.



THE ONLY 26 GAUGE
SAFETY CATHETER.
HELPING CLINICIANS
MEET THE NEEDS
OF THE SMALLEST
PATIENTS.



REFERENCES:

1. Delta Med data on file.
2. INS Standards of Practice, page S51, 26. B

For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Distributed by:
Smiths Medical International
1500 Eureka Park, Lower Pemberton
Ashford Kent, TN25 4BF
Tel: +44 (0)845 850 0445
www.smiths-medical.com

Manufactured by:
Delta Med SpA
Via G. Rossa, 20-46019
Viadana (Mn) - Italy
0375 785915

MPAUC-1685

**Rx
ONLY**

smiths medical

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. The product referenced is CE-marked. ©2017 Smiths Medical. All rights reserved. SS195458GB-102017