

**smiths medical**  
bringing technology to life



## DeltaVen® Closed System Catheter


Feel the difference,  
see the difference™



# FEEL THE DIFFERENCE, SEE THE DIFFERENCE™

Finding a catheter that can manage the rigor or power injection can be a challenge. DeltaVen® was designed to meet all your power injection needs.

A CATHETER DESIGNED TO TRAVEL WITH YOUR PATIENT THROUGHOUT THE HOSPITAL



- The only catheter rated to 330 PSI\* with flow rates up to 8 mL per second in a 20G catheter
- Helps eliminate the need for an additional catheter stick at radiology, CT or other areas
- Bonded extension set reduces the risk of contamination from repeated disconnections and reconnections

\*Please refer to Instructions For Use provided with the device for a full listing of the Indications and Precautions.

## PRESSURE RATED FOR 330 PSI IN GAUGES 24 THROUGH 16

IV Catheter Gauge	IV Catheter Length (mm)	Max Flow at 20°C (mL/sec)	Max Flow at 37°C (mL/sec)
24G	19	4	5
	19	5	6
22G	25	5	5
	25	7	8
20G	32	6	7
	45	6	7
18G	25	8	9
	32	7	9
16G	45	7	9
	32	>10	>10

Contrast media used: IOMERON 400 featuring 27,5 mPa.s. viscosity at 20°C and 12,6 mPa.s. at 37°C. Delta Med data on file.



For more information visit our website at [www.smiths-medical.com/deltaven](http://www.smiths-medical.com/deltaven)

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

**Distributed by:**  
Smiths Medical International  
1500 Eureka Park, Lower Pemberton  
Ashford Kent, TN25 4BF  
Tel: +44 (0)845 850 0445  
[www.smiths-medical.com](http://www.smiths-medical.com)

**Manufactured by:**  
**Delta Med SpA**  
Via G. Rossa, 20-46019  
Viadana (Mn) - Italy  
0375 785915

MPAUC-1686

**Rx  
ONLY**

smiths medical

Find your local contact information at: [www.smiths-medical.com/customer-support/contact-us](http://www.smiths-medical.com/customer-support/contact-us)

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. The product referenced is CE-marked. ©2017 Smiths Medical. All rights reserved. SS195459GB-102017