

smiths medical

level 1®



NORMOFLO® Irrigation Fluid Warmer

Minimize hypothermia and maintain comfort when delivering high volume irrigation.



THE NORMOFLO® IRRIGATION FLUID WARMER IS PERFECT FOR A VARIETY OF PROCEDURES INCLUDING:

- Laparoscopy
- Arthroscopy
- Hysteroscopy
- Wound debridement
- Urology
- Liposuction
- Lavage

The NORMOFLO® irrigation fluid warmer is designed to warm irrigating fluids during clinical scenarios. It is an important tool for infusing normothermic fluids and helping prevent Unintended Perioperative Hypothermia (UPH).

REDUCE BLEEDING, SPEED RECOVERY TIMES AND LOWER INFECTION RATES

With the NORMOFLO® irrigation fluid warmer, you can deliver large volumes of body temperature irrigation at optimum flow rates. This device helps you effectively maintain your patients' core temperature during the entire length of the procedure. Minimizing hypothermia has been shown to reduce bleeding, speed recovery and lower infection rates.^{1,2}

NO CHANGE IN SURGICAL TECHNIQUE

The NORMOFLO® irrigation fluid warmer delivers a clear, bubble-free field of vision, which assists surgeons to operate with confidence and with no change in surgical technique. It easily interfaces with standard surgical tools, hand pieces, scopes and specialized pumps to provide consistently safe performance that helps improve patient comfort and delivers better outcomes.

EASY TO SET UP AND USE

Model H-1100 is intended for gravity flow procedures and features an electrically operated IV pole that can automatically lift up to four 3-liter bags at the touch of a button.

ORDERING INFORMATION

Hardware		
Product Code	Description	Units per case
H-1100	NORMOFLO® irrigation fluid warmer includes electrically operated IV pole	1
Disposables		
IRI-600	Irrigating set, for all surgical irrigating gravity flow, or can be pressurized up to 300 mm/Hg.	10
IRI-600B	Irrigating set, for all surgical irrigating gravity flow, or can be pressurized up to 300 mm/Hg. Has spikes compatible with Baxter Healthcare's bags.	10
IR-700	Irrigating set, for all surgical irrigating, gravity flow, or can be pressurized up to 300 mm/Hg.	9
PL-6	Replacement lines for IRI-600	20

REFERENCES

1. Kurz A et al. Perioperative normothermia to reduce the incidence of surgical wound infection and shorten hospitalization. *The New England Journal of Medicine*. 1996;334(19): 1209-15.
2. Schmied H et al. Mild hypothermia increases blood loss and transfusion requirements during total hip arthroplasty. *Lancet*. 1996;347(8997):389-92.



The pressure system accommodates standard 1-liter bags in rigid chambers.

Lift heavy bags at the touch of a button



For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical ASD, Inc.
6000 Nathan Lane North
Minneapolis, MN, 55442, USA
Tel: 1-614-210-7300
Toll-Free USA 1-800-258-5361
www.smiths-medical.com

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Level 1, NORMOFLO and the Smiths Medical design mark are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2017 Smiths Medical. All rights reserved. TM194812EN-122017

MMSPCA-0078

**Rx
ONLY**

smiths medical