

Snuggle Warm® Tube Blanket



The convective warming blanket that surrounds the patient with warm air while providing unrestricted full body access

Key features of the Snuggle Warm® Tube Blanket

Flexible

- Ideal for use in the cath lab or any clinical setting where full body access is required
- Radiotranslucent material will not interfere with imaging procedures

Easy to Use

- Integrated straps secure blanket in place
- Compatible with the Level 1® EQUATOR® and Snuggle Warm® Convective Warmers


Enhanced Patient Satisfaction

- Promotes patient comfort by minimizing thermal discomfort and complaints of cold


SNUGGLE WARM® CONVECTIVE WARMING BLANKETS

Snuggle Warm® Blankets are...


- Made of soft comfortable fabric-like material — no hot plastic
- Strong and durable — resists tears, punctures and fluids
- Easy to manipulate — can be folded and repositioned for quick patient access
- Designed to conform to the patient's body with no billowing — no need to tie down or secure the blanket
- Latex free
- Radiotranslucent
- Available in multiple configurations for all clinical applications




Adult Full Body




Adult Lower Body




Adult Upper Body




Sterile Cardiac Blanket




Pediatric Full Body



Tube Blanket



Neonate to Small Child
and
Preemie to Neonate



Ordering Information

Ref. No.	Description	Units/Case
EQ-5000	EQUATOR® Convective Warming System	1
RC-5000	Rolling Cart	1

Snuggle Warm® Convective Warming Blankets

SW-2001	Adult Full Body Blanket	10
SW-2002	Pediatric Full Body Blanket	10
SW-2003	Adult Upper Body Blanket	10
SW-2004	Adult Lower Body Blanket	10
SW-2005	Neonate to Small Child Blanket	10
SW-2006	Preemie to Neonate Blanket	10
SW-2007	Tube Blanket	10
SW-2008	Sterile Cardiac Blanket	10

smiths

Smiths Medical
2231 Rutherford Road
Carlsbad, CA 92008
T: 800 848 1757
F: 800 621 2590
www.smiths-medical.com

Please refer to the Instructions for Use supplied with the product for detailed instructions, warnings, and cautions.

Level 1, EQUATOR, Controlled Convective Warming, Temperature Management Solutions, and Snuggle Warm are trademarks of the Smiths Medical family of companies. The symbol® indicates it is registered in the U.S. Patent and Trademark Office and certain other countries. The products described here are covered by one or more of the following U.S. Patent Nos.: 5,785,723; 6,143,020; 6,277,144; 6,440,157; 5,063,994; 5,097,898. Other patents pending

©2005 Smiths Medical family of companies. All rights reserved.

**Rx
ONLY**

TM-002 8/05