

# Warmer patients. Better outcomes.

Snuggle Warm® Convective Warming Lateral Access Blanket



TEMPERATURE MANAGEMENT


## Quality. Choice. Confidence. Now available in LATERAL ACCESS blankets.



- Proprietary design optimizes lower body procedures
- Durable material is resistant to tears, punctures, and fluids
- Fabric-like feel is soft and comfortable next to the patient's skin
- Precisely positioned air-holes and channel design to maximize heat delivery and keep patients warm

**SW-2014L/R Lateral Access Blankets are recommended for use in these procedures:**

- Hip Procedures
- Upper or Lower Leg Procedures



### Level 1® Ordering Information Snuggle Warm® Convective Warming Blankets (10/case):

**SW-2014L Left Lateral Access Blanket**  
**SW-2014R Right Lateral Access Blanket**



EQ-5000 Equator® Convective Warmer brings you patented temperature control ensuring delivery of desired therapy.

**Visit [www.smiths-medical.com/SnuggleWarm](http://www.smiths-medical.com/SnuggleWarm) for more information on convective warming.**

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

#### **Smiths Medical ASD, Inc.**

Norwell, MA 02061, USA  
Phone: 1-214-618-0218  
Toll-Free USA: 1-800-258-5361

[www.smiths-medical.com](http://www.smiths-medical.com)

**Smiths Medical is part of the global technology business Smiths Group plc.**

Product(s) described may not be licensed or available for sale outside of the United States.

Equator, Snuggle Warm, Level 1 and the Smiths Medical and Level 1 design mark are trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. © 2010 Smiths Medical family of companies. All right reserved. TM 194832 07/10

**smiths medical**