

# smiths medical

## level 1®



### Snuggle Warm® Pediatric Convective Warming Blankets

Warmth, comfort and security for your smallest patients. Snuggle Warm® convective warming blankets are designed to help maintain normothermic temperature in neonates, infants and children.



SWU-2002



SWU-2009



SWU-2011

## SNUGGLES THE PATIENT WHEN INFLATED

Snuggle Warm® convective warming blankets help provide the warmth, comfort and security your smallest patients need when they are most vulnerable. Made of soft fabric-like material that is comfortable next to the skin, the blankets conform to patient anatomy to maximize heat retention. The unique durable construction also resists tears and punctures. On pediatric underbody blankets, absorbent pads minimize pooling of fluid, and an uninflated headrest accommodates head pads.

## SNUGGLE WARM® PEDIATRIC BLANKET SPECIFICATIONS

Product Number	Name	Size (uninflated)	Weight	Case Quantity	Hose Ports
SWU-2002	Pediatric full body blanket	102 cm W x 146 cm L (40" W x 57.5" L)	160 g (5.6 oz)	10	1
SWU-2009	Pediatric underbody blanket	66 cm W x 94 cm L (26" W x 37" L)	120 g (4.0 oz)	10	2
SWU-2011	Large pediatric underbody blanket	102 cm W x 136 cm L (40" W x 53" L)	200 g (7.2 oz)	10	2



Blanket model SWU-2011

ASK ME ABOUT  
LEVEL 1®  
TEMPERATURE  
MONITORING PROBES!

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

MPAUC-1179

**Smiths Medical ASD, Inc.**  
St. Paul, MN 55112, USA  
Phone: 1-614-210-7300  
Toll-Free USA: 1-800-258-5361  
[www.smiths-medical.com](http://www.smiths-medical.com)

**Rx  
ONLY**  
smiths medical

Find your local contact information at: [www.smiths-medical.com/customer-support](http://www.smiths-medical.com/customer-support)

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Snuggle Warm, Level 1, HOTLINE, Neo-Therm and the Level 1 and Smiths Medical design marks are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. ©2016 Smiths Medical. All rights reserved. TM194867EN-062019