

Portex® ULTRApertc® Percutaneous Dilation Tracheostomy Kits with Single Stage Curved Dilator



Full Component or Basic Kits—
with or without Portex® Blue
Line ULTRA® or Portex® Blue
Line ULTRA® Suctionaid®
Tracheostomy Tubes.



Pre-lubricated
ergonomic Single Stage
curved dilator.



Introducer with
handle included
in kits with tubes.



Introducer compatible
with Portex® Blue Line
ULTRA® and Portex® Blue
Line ULTRA® Suctionaid®
Tracheostomy Tubes.

The Portex® ULTRAperc® Percutaneous Dilation Tracheostomy Kits allow the percutaneous insertion of a tracheostomy tube using the widely accepted Seldinger technique by means of a Single Stage curved dilator. Kits are available with or without the Portex® Blue Line ULTRA® or Portex® Blue Line ULTRA® Suctionaid® Tracheostomy Tubes.

These kits feature a Single Stage, pre-lubricated curved dilator with ergonomic handle that facilitates insertion. The specially designed introducer can also be used to percutaneously insert the Portex® Blue Line ULTRA® or Portex® Blue Line ULTRA® Suctionaid® Tracheostomy Tubes.

ULTRAperc® Full Kit Components

- | | | | |
|---|---|--------------------------|------------------------------------|
| (2) Lidocaine HCl (1.5%) and Epinephrine (1:200,000) 5 ml | (1) Scalpel with #15 Blade | (8) Gauze Sponge | (2) Povidone-Iodine Swabstick, 3's |
| (1) Guide Wire Assembly | (1) Pre Dilator | (2) Point-Lok® Device | (1) Suture Kit |
| (1) 14 G Needle and Cannula Assembly | (1) Single Stage Dilator | (1) 22 G x 1½" Needle | (1) Curved Forceps, 5" |
| (1) Guiding Catheter | (3) Plastic Syringes (10 ml, Luer Slip) | (1) Fenestrated Drape | |
| | | (1) Split Trach Dressing | |
| | | (1) Swivel Adapter | |

| REFERENCE CODE | DESCRIPTION | TUBE I.D. | TUBE O.D. | TUBE LENGTH | 2 UNITS PER CASE |
|----------------|--|-----------|-----------|-------------|------------------|
| 100/596/070 | 7.0 mm Portex® Blue Line ULTRA® Suctionaid® Kit | 7.0 Size | 10.5 mm | 70.0 mm | |
| 100/596/080 | 8.0 mm Portex® Blue Line ULTRA® Suctionaid® Kit | 8.0 Size | 11.9 mm | 75.5 mm | |
| 100/596/090 | 9.0 mm Portex® Blue Line ULTRA® Suctionaid® Kit | 9.0 Size | 13.3 mm | 81.0 mm | |
| 100/595/070 | 7.0 mm Portex® Blue Line ULTRA® Kit | 7.0 Size | 10.5 mm | 70.0 mm | |
| 100/595/080 | 8.0 mm Portex® Blue Line ULTRA® Kit | 8.0 Size | 11.9 mm | 75.5 mm | |
| 100/595/090 | 9.0 mm Portex® Blue Line ULTRA® Kit | 9.0 Size | 13.3 mm | 81.0 mm | |
| 100/573/000 | Kit with Single Stage Dilator and Soft Introducers for 7, 8 and 9 mm Tubes Only; No Tracheostomy Tube, Trach Wedge or Trach Tube Holder with Brush Included with this Kit. | | | | |

ULTRAperc® Basic Kit Components

- | | | | |
|--------------------------------------|---|-------------------------|----------------------|
| (1) 14 G Needle and Cannula Assembly | (1) Plastic Syringes (10 ml, Luer Slip) | (1) Pre Dilator | (1) Guiding Catheter |
| (1) Scalpel with #15 Blade | (1) Single Stage Dilator | (1) Guide Wire Assembly | (8) Gauze Sponge |

| REFERENCE CODE | DESCRIPTION | TUBE I.D. | TUBE O.D. | TUBE LENGTH | 1 UNIT PER CASE |
|----------------|--|-----------|-----------|-------------|-----------------|
| 100/591/070 | 7.0 mm Portex® Blue Line ULTRA® Kit | 7.0 Size | 10.5 mm | 70.0 mm | |
| 100/591/080 | 8.0 mm Portex® Blue Line ULTRA® Kit | 8.0 Size | 11.9 mm | 75.5 mm | |
| 100/591/090 | 9.0 mm Portex® Blue Line ULTRA® Kit | 9.0 Size | 13.3 mm | 81.0 mm | |
| 100/572/000 | Kit with Single Stage Dilator and Soft Introducers for 7, 8 and 9 mm Tubes Only; No Tracheostomy Tube, Trach Wedge or Trach Tube Holder with Brush Included with this Kit. | | | | |

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

Smiths Medical ASD, Inc.
 160 Weymouth Street, Rockland, MA 02370
 Phone: (800) 848-1757 Fax: (800) 621-2590
 www.smiths-medical.com
Smiths Medical, part of the global technology business Smiths Group.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician. Please refer to the instructions for use supplied with the product for detailed instructions, warnings, and cautions. Smiths Medical design mark, Portex, Blue Line ULTRA, and ULTRAperc are registered trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.
 ©2008 Smiths Medical family of companies. All rights reserved. RE453 - TR194420EN 10.2008

