

# For Optimal Fit

## Adjustable for Better Fit


In the past, doctors had no choice but to use ET tubes for patients with thick, bull necks. Now, there is a tube with an adjustable flange that will allow clinicians to achieve optimal fit in patients that require extra length (horizontal and/or vertical) dimensions.

## The Silicone Story

The unique biocompatibility characteristics of silicone offer clinicians and patients a more tissue friendly option.

Additional benefits include:

- Durable – less affected by time, light, temperature
- Remains naturally soft and flexible
- Hydrophobic-resistant to water absorption
- Latex-free



### ORDERING INFORMATION

#### TIGHT TO SHAFT (TTS™) ADJUSTABLE NECK FLANGE HYPERFLEX™ TRACHEOSTOMY TUBES

Reference Code	Inside Diameter	Outside Diameter	Tube Length	Cuff Resting Diameter
67HA60	6.0 mm	9.2 mm	110.0 mm	11.0 mm
67HA70	7.0 mm	10.6 mm	120.0 mm	12.0 mm
67HA80	8.0 mm	11.7 mm	130.0 mm	13.0 mm
67HA90	9.0 mm	12.9 mm	140.0 mm	14.5 mm

Supplied sterile; individually packaged with clip-in, curved obturator with hollow tip, twill trach tie and decannulation cap.

#### MID-RANGE AIRE-CUF® ADJUSTABLE NECK FLANGE HYPERFLEX™ TRACHEOSTOMY TUBES

Reference Code	Inside Diameter	Outside Diameter	Tube Length	Cuff Resting Diameter
75HA60	6.0 mm	9.2 mm	110.0 mm	18.0 mm
75HA70	7.0 mm	10.6 mm	120.0 mm	21.0 mm
75HA80	8.0 mm	11.7 mm	130.0 mm	25.0 mm
75HA90	9.0 mm	12.9 mm	140.0 mm	28.0 mm


Supplied sterile; individually packaged with clip-in, curved obturator with hollow tip and twill trach tie.

# Versatile

For use in the hospital and home

- Fixed Flange
- Less defined curvature allows for fitting multiple tracheal anatomies
- Can be resterilized up to 10 times
- Longer length tube than the most commonly used extra length tube

Flexible, wire-reinforced shaft can adapt to unique anatomical needs.



## ORDERING INFORMATION

### UNCUFFED FIXED FLANGE HYPERFLEX™ EXTRA LENGTH TRACHEOSTOMY TUBES

Reference Code	Inside Diameter	Outside Diameter	Tube Length	Cuff Resting Diameter
60AFHXL60	6.0 mm	9.2 mm	100.0 mm	NA
60AFHXL70	7.0 mm	10.6 mm	110.0 mm	NA
60AFHXL80	8.0 mm	11.7 mm	120.0 mm	NA
60AFHXL90	9.0 mm	12.9 mm	130.0 mm	NA

Supplied sterile; individually packaged with clip-in, curved obturator with hollow tip, twill trach tie and decannulation cap.

### CUFFED FIXED FLANGE HYPERFLEX™ EXTRA LENGTH TRACHEOSTOMY TUBES

Reference Code	Inside Diameter	Outside Diameter	Tube Length	Cuff Resting Diameter
75FHXL60	6.0 mm	9.2 mm	100.0 mm	18.0 mm
75FHXL70	7.0 mm	10.6 mm	110.0 mm	21.0 mm
75FHXL80	8.0 mm	11.7 mm	120.0 mm	25.0 mm
75FHXL90	9.0 mm	12.9 mm	130.0 mm	28.0 mm

Supplied sterile; individually packaged with clip-in, curved obturator with hollow tip and twill trach tie.

### TIGHT TO SHAFT (TTS™) CUFFED EXTRA LENGTH FIXED FLANGE HYPERFLEX™ TRACHEOSTOMY TUBES

Reference Code	Inside Diameter	Outside Diameter	Tube Length	Cuff Resting Diameter
67FHXL60	6.0 mm	9.2 mm	100.0 mm	11.0 mm
67FHXL70	7.0 mm	10.6 mm	110.0 mm	12.0 mm
67FHXL80	8.0 mm	11.7 mm	120.0 mm	13.0 mm
67FHXL90	9.0 mm	12.9 mm	130.0 mm	14.5 mm

Supplied sterile; individually packaged with clip-in, curved obturator with hollow tip, twill trach tie and decannulation cap.

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

**Smiths Medical ASD, Inc.**  
Norwell, MA 02061, USA  
Phone: 1-214-618-0218  
Toll-Free USA: 1-800-258-5361

**Smiths Medical Canada Ltd.**  
Markham, Ontario, Canada, L3R 4Y8  
Phone: 905-477-2000  
Toll-Free: 1-800-387-4346

Customer and Clinical Services  
U.S. 1-800-258-5361  
Canada 1-800-387-4346  
www.smiths-medical.com

www.smiths-medical.com

Smiths Medical is part of the global technology business Smiths Group plc.

Product(s) described may not be licensed or available for sale outside of the United States.

Portex and Smiths Medical design marks, TTS, Hyperflex and Aire-Cuf are trademarks of the Smiths Medical family of companies.

The symbol © indicates that the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.

All other names and marks mentioned are trade names, trademarks, and service marks of their respective owners.

©2010 Smiths Medical family of companies. All rights reserved. TR194436 EN 03.2010

**smiths medical**