

## SurgiVet® Hand Held Pulse Oximeter



### VETERINARY

- Low perfusion performance
- Motion tolerance
- Pulse amplitude Index (PI) bar-graph indicator to help identify and maintain a good test site
- Quick selection of default alarm limit values by patient type: Cat, Dog, Horse, Factory and Custom
- Large, bright LED display
- Memory storage for up to 99 patients
- Audible and visual alarms
- Multiple power choices
- Three-year warranty

**SurgiVet**<sup>TM</sup>

## SurgiVet® V1030 Hand Held Pulse Oximeters

Built on a foundation of proprietary technology, the new SurgiVet® V1030 Hand Held Pulse Oximeter provides fast, reliable SpO<sub>2</sub> pulse rate and pulse signal strength measurements across a full range of patients and helps maintain accurate SpO<sub>2</sub> readings during periods of patient motion and when monitoring patients with low perfusion.



V1030 pulse oximeter monitor system includes:

- 5 ft. oximetry extension cable
- Reflectance oximetry sensor
- Universal “Y” oximetry sensor
- Four (4) “AA” alkaline batteries (17 hrs of continuous use)
- Blue Protective Glove
- A CD that contains the Operation Manual, Service Manual, and Quick Reference Guide. A printed Quick Reference Guide is also included.



WW1026SYS

- Download real-time data, graphical and numerical trends
- Modular style allows for attachment to Oximeter Docking Station
- Includes 5 Pack of printer paper



WW1025SYS

- Converts handheld device to tabletop unit
- Allows for convenient storage of device
- Recharges device and additional backup Lithium Ion Battery
- Modular style allows for attachment of optional printer

Catalog #	Description
<b>Power Supplies</b>	
WW1090	Custom Lithium Ion Rechargeable Battery (30 hrs continuous use; 3hr recharge <sup>1</sup> )
WW1095	Power Supply Universal with US/Japan Adapter 30 WATT, 100-240V, 50Hz
WW1095OUS	Power Supply Universal with International Blades 30 WATT, 100-240V, 50Hz
WW1096SYS	AC Power and Rechargeable Lithium Ion Battery System
WW1096OUS	AC Power and Rechargeable Lithium Ion Battery System (OUS*)
<b>Mounting Brackets</b>	
WW1098	Universal Mounting Bracket
WW1098C	Universal Mounting Bracket with IV Pole Clamp
<b>Attachable Printer</b>	
WW1026SYS	Attachable Thermal Printer (attaches to the docking station)
WW1027	Printer Paper, 5 Pack
<b>Docking Stations</b>	
WW1025SYS	Docking Station
WW1025OUS	Docking Station (OUS*)
WW1028SYS	Docking Station and Printer System
WW1028OUS	Docking Station and Printer System (OUS*)

<sup>1</sup>Approximation based on typical usage. \*OUS = Outside United States

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE. PLEASE READ THE INSTRUCTIONS FOR USE SUPPLIED WITH THE PRODUCT FOR DETAILED INSTRUCTIONS, WARNINGS AND CAUTIONS.

**Smiths Medical ASD, Inc.**  
 St. Paul, MN 55112, USA  
 Phone: 1-214-618-0218  
 Toll-Free USA: 1-800-258-5361

**Customer and Clinical Services**  
 USA: 1-214-618-0218  
 Toll-Free USA: 1-800-258-5361

**Customer Service**  
 Canada 1-905-477-2000  
 Toll-Free Canada: 1-800-387-4346  
[www.surgivet.com](http://www.surgivet.com)  
[www.smiths-medical.com](http://www.smiths-medical.com)

**Smiths Medical is part of the global technology business Smiths Group plc.**

SurgiVet and the Smiths Medical design marks are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.  
 ©2010 Smiths Medical. All rights reserved. VE195816 01.2011