

smiths medical  
medex™



## Clear-Cuff® Pressure Infusor Bags

Pressure Monitoring




# Clear-Cuff® Pressure Infusor Bags

## Durable

Durable, easy-to-clean, polyurethane construction allows the Clear-Cuff® pressure infusor bags to be re-used if desired.

## Circumferential Bladder Construction<sup>1</sup>

The Clear-Cuff® pressure infusor bags provide a consistent flow rate for effective fluid management. Its circumferential bladder construction allows for an optimal squeeze of the entire fluid bag, enabling more consistent flow rates.



## Beneficial Features

The Clear-Cuff® pressure infusor bags, made of thick, durable polyurethane, have an inherent transparent quality, providing easy viewing of the fluid bag and the fluid level. Additionally, a special built-in hook allows the fluid bag to be replaced without removing it from the IV pole. The Clear-Cuff® pressure infusor bags also feature an extra-large squeeze bulb and pressure gauge for fast and easy inflation.



Clear-Cuff® Pressure Infusor Bags (D Fusors)		
Reference Code	Description	Units Per Case
MX4705P1CZ	500 mL Pressure Infusor	10
MX4710P1CZ	1000 mL Pressure Infusor	10

## Reference

1. Data on file

For more information visit our website at [www.smiths-medical.com](http://www.smiths-medical.com)

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

**Smiths Medical India Private Limited**  
 508/509, Western Edge-2, W.E. Highway,  
 Borivali (E), Mumbai - 66, India.  
 T: + 99 22 28705210 F: +91 22 28705211  
 E: [csinda@smiths-medical.com](mailto:csinda@smiths-medical.com)  
[www.smiths-medical.com](http://www.smiths-medical.com)

Find your local contact information at: [www.smiths-medical.com/customer-support](http://www.smiths-medical.com/customer-support)

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Clear-Cuff, Medex and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2019 Smiths Medical. All rights reserved.  
 VA0020.GB.IN.Rev.A.0919

**Rx  
ONLY**

smiths medical