

smiths medical
deltec®



SAFESHEATH® ULTRA LITE

Valved Peel-Away Introducers




SAFESHEATH® ULTRA LITE VALVED PEEL-AWAY INTRODUCERS

Designed for lower risk of blood loss and air embolisms

Smiths Medical now offers a new valved peel-away introducer designed for lower risk of blood loss and air embolisms. Available in a wide variety of sizes and for use with PORT-A-CATH® implantable access systems.

- Low insertion force valve membrane allows smooth insertion and removal of port catheters
- Extruded scoreline peel-away sheath for clean, even peeling during removal
- Latex-free



INFORMATION

SafeSheath® introducers are available in Smiths Medical PORT-A-CATH® system trays and also as separate sets.

PORT-A-CATH® II Systems with SafeSheath® Valved Peel-Away Introducers						
Product	Introducer Size	Catheter Material	O.D. (mm)	Length (cm)	Connector	Tray Reorder No.
Single-lumen portal	8.5 Fr	Polyurethane	2.6	76	ULTRA-LOCK®	21-4066-24
Single-lumen portal	6 Fr		1.9	76	ULTRA-LOCK®	21-4067-24
Single-lumen Low Profile™ portal	8.5 Fr		2.6	76	ULTRA-LOCK®	21-4068-24
Single-lumen Low Profile™ portal	7 Fr		1.9	76	ULTRA-LOCK®	21-4069-24

PORT-A-CATH® POWER P.A.C. Systems with SafeSheath® Valved Peel-Away Introducers						
Product	Introducer Size	Catheter Material	O.D. (mm)	Length (cm)	Connector	Tray Reorder No.
Single-lumen portal	8.5 Fr	Polyurethane	2.6	76	ULTRA-LOCK®	21-4474-24
Single-lumen portal	6 Fr		1.9	76	ULTRA-LOCK®	21-4476-24
Single-lumen Low Profile™ portal	8.5 Fr		2.6	76	ULTRA-LOCK®	21-4475-24
Single-lumen Low Profile™ portal	7 Fr		1.9	76	ULTRA-LOCK®	21-4477-24
Dual-lumen Low Profile™ portal	7 Fr		2.2	76	WING-LOCK®	21-8469-24
Dual-lumen Low Profile™ portal	10 Fr		3.2	76	WING-LOCK®	21-8470-24

SafeSheath® Valved Peel-Away Introducer Sets	
Introducer Sets (5/box)	
French Size	Reorder #
6 Fr	21-2476-24
7 Fr	21-2477-24
8.5 Fr	21-2478-24
9 Fr	21-2479-24
10 Fr	21-2480-24
11 Fr	21-2481-24

Introducer sheath length – 13 cm
Dilator sheath length – 18 cm

For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES.

Smiths Medical ASD, Inc.
6000 Nathan Lane North
Minneapolis, MN, 55442, USA
Tel: 1-614-210-7300
Toll-Free USA 1-800-258-5361
www.smiths-medical.com

MKEECA-0191

**Rx
ONLY**

smiths medical

Smiths Medical is part of the global technology business Smiths Group plc. Product(s) described may not be licensed or available for sale in all countries. PORT-A-CATH, ULTRA-LOCK, WING-LOCK, Low Profile, and the Smiths Medical and Deltec design marks are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. Product referenced is CE-Marked ©2018 Smiths Medical. All rights reserved. VA192200EN-052018