


smiths medical deltec®



GRIPPER PLUS® POWER P.A.C. Safety Huber Needles

Proven performance you expect,
now with power capabilities



GRIPPER PLUS® POWER P.A.C. needles provide the flexibility of use for any therapy need, including power injection of contrast media. GRIPPER PLUS® POWER P.A.C. needles allow for higher pressure and flow rates for power injections with PORT-A-CATH® POWER P.A.C. systems.


GRIPPER PLUS® POWER P.A.C. NEEDLES OFFER:

- Stable, cushioned platform
- Strong, latex-free plastic tubing
- Protected access site
- Needlestick protection
- Luer-lock connector
- 5 ml/sec flow rate

GRIPPER PLUS® POWER P.A.C. Safety Huber Needles

The GRIPPER PLUS® POWER P.A.C. needle's unique de-accessing method lifts the needle straight out with a single, smooth motion. The needle then "clicks" into a locked position, which protects you from needlesticks.

I.V. catheters may vary greatly with regard to materials, performance, and safety features. This guide is designed primarily for cross-referencing I.V. catheter gauges and lengths. Please contact your local Smiths Medical representative for assistance in cross-referencing to the I.V. catheter that best fits your clinical needs.



ORDERING INFORMATION

GRIPPER PLUS® POWER P.A.C. Safety Huber Needles				
Gauge	Length	Maximum Flow Rate	Reorder Numbers NAC	Reorder Numbers Non Y-site
19	3/4 inch	5ml/sec	21-3468-24	21-3368-24
19	1 inch	5ml/sec	21-3469-24	21-3364-24
19	1 1/4 inch	5ml/sec	21-3470-24	21-3365-24
20	5/8 inch	5ml/sec	21-3464-24	21-3369-24
20	3/4 inch	5ml/sec	21-3465-24	21-3367-24
20	1 inch	5ml/sec	21-3466-24	21-3362-24
20	1 1/4 inch	5ml/sec	21-3467-24	21-3363-24

For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical Australasia Pty Ltd
Suite 201, Level 2
10 Norbrik Drive
Bella Vista NSW 2153, Australia
Tollfree- Australia 1800 654 949
www.smiths-medical.com

Distributed by:
USL Medical
494 Rodebank Road
Avondale, 1026 New Zealand
Freecall: 0800 658 814

MMSPCA-0155

smiths medical

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. GRIPPER PLUS, PORT-A-CATH, P.A.S. PORT, Deltec and the Smiths Medical design are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. ©2018 Smiths Medical. All rights reserved. VA192500GBAUS-022018