



100/100/055

AIRCARE™ PVC Endotracheal Tubes

- Medical grade PVC is transparent and softens at body temperature to provide patient comfort
- Radiopaque blue line aids in confirmation of correct tube placement via x-ray
- Smooth tip helps to minimize trauma during oral intubation
- 15 mm connector included with each tube
- Murphy eye to minimize potential occlusion of the tube
- 3.0 mm to 10.0 mm including half sizes
- An effective alternative for short term intubation
- Used for airway management
- Aid in mechanical ventilation
- A low cost alternative to help prevent disease transmission
- Latex free
- Sterile
- Single Use*

*The following changes have recently been added to the AAHA Standards of Accreditation:

ANESTHESIA:

AN30.1. Endotracheal tubes are single-use disposable, autoclaved or gas sterilized between each use.

Catalog #	Description	Size (I.D.)	Size (O.D.)	Length	Units/Box
100/100/030	3.0 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	3.0 mm	4.5 mm	16.5 cm	10
100/100/035	3.5 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	3.5 mm	5.3 mm	18.5 cm	10
100/100/040	4.0 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	4.0 mm	6.0 mm	21.0 cm	10
100/100/045	4.5 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	4.5 mm	6.5 mm	22.5 cm	10
100/100/050	5.0 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	5.0 mm	6.9 mm	24.5 cm	10
100/100/055	5.5 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	5.5 mm	7.5 mm	27.5 cm	10
100/100/060	6.0 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	6.0 mm	8.2 mm	28.5 cm	10
100/100/065	6.5 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	6.5 mm	8.7 mm	29.5 cm	10
100/100/070	7.0 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	7.0 mm	9.3 mm	31.0 cm	10
100/100/075	7.5 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	7.5 mm	10.0 mm	31.5 cm	10
100/100/080	8.0 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	8.0 mm	10.7 mm	33.0 cm	10
100/100/085	8.5 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	8.5 mm	11.3 mm	33.0 cm	10
100/100/090	9.0 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	9.0 mm	12.0 mm	33.0 cm	10
100/100/010	10.0 mm AIRCARE™ Cuffed Endotracheal Tube (Sterile)	10.0 mm	13.3 mm	33.0 cm	10

Smiths Medical ASD, Inc.

Norwell, MA 02061, USA
Phone: 1-214-618-0218
Toll-Free USA: 1-800-258-5361
www.smiths-medical.com
www.surgivet.com

Smiths Medical, part of the global technology business Smiths Group

Product(s) described may not be licensed or available for sale in Canada or other countries outside of the United States.

AIRCARE and the Smiths Medical and SurgiVet design marks are trademarks of the Smiths Medical family of companies. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners. © 2009 Smiths Medical family of companies. All rights reserved. AIRCARE Rev.01 09/09

