

SurgiVet® Hand Held Pulse Oximeter



VETERINARY

- Low perfusion performance
- Motion tolerance
- Pulse amplitude Index (PI) bar-graph indicator to help identify and maintain a good test site
- Quick selection of default alarm limit values by patient type: Cat, Dog, Horse, Factory and Custom
- Large, bright LED display
- Memory storage for up to 99 patients
- Audible and visual alarms
- Multiple power choices
- Three-year warranty

SurgiVetTM

SurgiVet® V1030 Hand Held Pulse Oximeters

Built on a foundation of proprietary technology, the new SurgiVet® V1030 Hand Held Pulse Oximeter provides fast, reliable SpO₂ pulse rate and pulse signal strength measurements across a full range of patients and helps maintain accurate SpO₂ readings during periods of patient motion and when monitoring patients with low perfusion.



V1030 pulse oximeter monitor system includes:

- 5 ft. oximetry extension cable
- Reflectance oximetry sensor
- Universal “Y” oximetry sensor
- Four (4) “AA” alkaline batteries (17 hrs of continuous use)
- Blue Protective Glove
- A CD that contains the Operation Manual, Service Manual, and Quick Reference Guide. A printed Quick Reference Guide is also included.



WW1026SYS

- Download real-time data, graphical and numerical trends
- Modular style allows for attachment to Oximeter Docking Station
- Includes 5 Pack of printer paper



WW1025SYS

- Converts handheld device to tabletop unit
- Allows for convenient storage of device
- Recharges device and additional backup Lithium Ion Battery
- Modular style allows for attachment of optional printer

Catalog #	Description
Power Supplies	
WW1090	Custom Lithium Ion Rechargeable Battery (30 hrs continuous use; 3hr recharge ¹)
WW1095	Power Supply Universal with US/Japan Adapter 30 WATT, 100-240V, 50Hz
WW1095OUS	Power Supply Universal with International Blades 30 WATT, 100-240V, 50Hz
WW1096SYS	AC Power and Rechargeable Lithium Ion Battery System
WW1096OUS	AC Power and Rechargeable Lithium Ion Battery System (OUS*)
Mounting Brackets	
WW1098	Universal Mounting Bracket
WW1098C	Universal Mounting Bracket with IV Pole Clamp
Attachable Printer	
WW1026SYS	Attachable Thermal Printer (attaches to the docking station)
WW1027	Printer Paper, 5 Pack
Docking Stations	
WW1025SYS	Docking Station
WW1025OUS	Docking Station (OUS*)
WW1028SYS	Docking Station and Printer System
WW1028OUS	Docking Station and Printer System (OUS*)

¹Approximation based on typical usage. *OUS = Outside United States

SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE. PLEASE READ THE INSTRUCTIONS FOR USE SUPPLIED WITH THE PRODUCT FOR DETAILED INSTRUCTIONS, WARNINGS AND CAUTIONS.

Smiths Medical ASD, Inc.

St. Paul, MN 55112, USA
Phone: 1-214-618-0218
Toll-Free USA: 1-800-258-5361

Customer and Clinical Services

USA: 1-214-618-0218
Toll-Free USA: 1-800-258-5361

Customer Service

Canada 1-905-477-2000
Toll-Free Canada: 1-800-387-4346

www.surgivet.com
www.smiths-medical.com

Smiths Medical is part of the global technology business Smiths Group plc.

SurgiVet and the Smiths Medical design marks are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.
©2011 Smiths Medical. All rights reserved. VE195816 01.2011