

smiths medical
animal health



Animal Consumables

ORDERING INFORMATION



Dog Urinary Catheter Luers

Product Code	Size (Fg)	Units per Case
303/800/060	6	50
303/800/080	8	50
303/800/100	10	50

Each catheter is designed with quality grade nylon. This allows for flexibility without rigidity for ease of insertion with minimal irritation.

Dog Oesophageal Tube

Product Code	Size (Fg)	Units per Case
403/810/230	23	40
403/810/350	35	40

Each tube is crafted with flexible siliconized Polyvinyl Chloride (PVC) tubing. Each tube comes with a rounded locking rod for ease of insertion with minimal irritation.



Cat Urinary Catheter Luers

Product Code	Size (Fg)	Units per Case
303/820/030	3	50
303/820/040	4	50

Each catheter is designed with quality grade nylon. This allows for flexibility without rigidity for ease of insertion with minimal irritation.



Cat Urinary Jackson Catheter

Product Code	Size (Fg)	Units per Case
303/825/030	3	50
303/825/040	4	50


Each catheter is designed with quality grade nylon. The Luer mount stylet provides rigidity without irritation for ease of insertion.



Horse Catheters

Product Code	Size (Fg)	Units per Case
303/840/200	20	100
303/840/240	24	100

Each catheter is designed with quality transparent Polyvinyl Chloride (PVC) tubing. The tip has a rounded locking rod for ease of insertion with minimal irritation.



Mare Catheter

Product Code	Size (Fg)	Units per Case
303/850/260	26	100

Each catheter is designed with quality transparent Polyvinyl Chloride (PVC) tubing. The tip has a rounded locking rod for ease of insertion with minimal irritation.

Tiemann Bitch Catheters

Product Code	Size	Units per Case
309/120/080	Special	10
309/120/100	10	10
309/120/120	12	10

Each catheter is designed with quality transparent Polyvinyl Chloride (PVC) tubing. The tip is curved for ease of insertion with minimal irritation.



Horse Stomach Tubes

Product Code	Size	Units per Case
403/800/100	Foal	100
403/800/200	Small	60
403/800/300	Medium	30
403/800/400	Large	20

Each tube is crafted with flexible siliconized Polyvinyl Chloride (PVC) tubing. Our tubes are designed for a variety of climates. Each tube comes with a funnel end to make the attachment to a continuous flow drench pump as quick and easy as possible.



Portex® Irrigating Cannula

Product Code	Color	Units per Case
413/820/010	Blue	50
413/820/020	Pink	50

Each catheter is designed with quality grade nylon. This allows for flexibility without rigidity for ease of insertion with minimal irritation.



For more information visit www.smiths-medical.com/surgivet

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical ASD, Inc.
6000 Nathan Lane North
Minneapolis, MN 55442, USA
Phone: 1-614-210-7300
Toll-Free USA 1-800-258-5361
www.smiths-medical.com

MHYTCA-1280

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us

smiths medical

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Medfusion, Surgivet and the Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2019 Smiths Medical. All rights reserved. VE195886EN-012019