

smiths medical
portex[®]

General Anaesthesia Airway

Breathing Systems Catalogue



Breathing
Circuits



Face
Masks



Circuit
Accessories





Contents

Expandable Length Dual-Limb Circuits	4
Fixed Length Dual-Limb Circuits.....	7
Fixed Length Single-Limb Circuits	10
Custom Circuits	11
VALUE PACK Custom Kits.....	12
Premium Soft Plus™ Face Masks with Inflation Valve	13
Breathing Bags	14
Bacterial / Viral Filter	15
Heat and Moisture Exchangers (HMEs).....	16
Filtered Heat and Moisture Exchangers (HMEFs).....	17
Accessories	18
Index.....	19

Expandable Length Dual-Limb Circuits

ORDERING INFORMATION

2x Adult 42"-120" (107 cm-305 cm) Expandable Length Dual-Limb Breathing Circuit Elbow: GSP Wye: Straight						
Item Number	Temp Port	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
450950-NL		Breathing	GSL without filter	3.0L (NL)	5245	12
450951-NL	Yes	Breathing (2), HME	GSL without filter	3.0L (NL)	5245	12
450952-NL		Breathing (2)		3.0L (NL)		12
451040-NL		HEPA		3.0L (NL)		12

2x Adult 42"-120" (107 cm-305 cm) Expandable Length Dual-Limb Breathing Circuit Elbow: GSP Wye: Straight				
Item Number	Tubing	Filter(s)	Breathing Bag	Units per Case
370900-NL	2x144" (366 cm)	Breathing (2)	3.0L (NL)	10

2x Adult 42"-108" (107 cm-274 cm) Expandable Length Dual-Limb Breathing Circuit Elbow: GSP Wye: Straight						
Item Number	Temp Port	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
450953-NL		Breathing	GSL no filter	3.0L (NL)	5245	20
450954-NL		Breathing		3.0L (NL)		20
451041-NL		HEPA	GSL no filter	3.0L (NL)		12
450955-NL		Breathing (2)		3.0L (NL)		20
450956-NL		Breathing, HME	GSL no filter	3.0L (NL)	5245	12

2x Pediatric 60"-150" (153-381 cm) MRI Length Dual-Limb Breathing Circuit Elbow: GSP Wye: Straight						
Item Number	Temp Port	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
371001-NL			GSL no filter	2.0L (NL)		10
371002			GSL no filter			10

2x Adult 120" (304 cm) Adult MRI Length Dual-Limb Breathing Circuit		
Item Number	Tubing	Units per Case
487000	2x 20"-40"-60" (51-102-152 cm)	20



ORDERING INFORMATION

2x Adult 33"-87" (84-221 cm) Expandable Length Dual-Limb Breathing Circuit						
Elbow: GSP Wye: Straight						
Item Number	Temp Port	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
450960-NL		Breathing	GSL without filter	3.0L (NL)	5245	20
450961-NL		Breathing (2)		3.0L (NL)		20
450962-NL		Breathing	GSL without filter	3.0L (NL)	005445	20
450963-NL		Breathing	GSL without filter	3.0L (NL)		20
451050-NL		HEPA		3.0L (NL)		20
450964-NL		Breathing (2)	GSL without filter	3.0L (NL)	5245	20
451051-NL		HME	GSL without filter	3.0L (NL)	5245	20
450965-NL	Yes	HME, Breathing		3.0L (NL)	5245	20
450904		BV		3.0L (L)		20
650904		BV		3.0L (NL)		20
450904-SG		BV		3.0L (L) Suregrip™ Bag		20
450907		BV		2.0L (L)		20
450907-NL		BV		2.0L (NL)		20
450906	Yes	BV		3.0L (L)		20
450906-NL	Yes	BV		3.0L (NL)		20
451004				3.0L (L)		20
451004-NL				3.0L (NL)		20
Elbow: GSP WYE: Bifurcated with port						
454204		BV		3.0L (L)		20
Elbow: GSP WYE: Parallel with Port						
455477D	Yes	BV	GSL without filter	3.0L (L)		20
455417D	Yes	BV	GSL with filter	3.0L (L)		20
Elbow: Std. WYE: Bifurcated with port						
455300				3.0L (L)		20
Elbow: - WYE: Swivel						
454000				3.0L (L)		20
453600		BV		3.0L (L)		20
4509901D		Breathing (2)	GSL without filter	3.0L (L)		20
Elbow: GSP WYE: Straight with port						
453600-NL		BV		3.0L (NL)		20



ORDERING INFORMATION

2x Pediatric 40"-100" (102-254 cm) Expandable Length Dual-Limb Breathing Circuit Elbow: GSP Wye: Straight					
Item Number	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
490900-NL	600 HME, Breathing	GSL without filter	1.0L (NL)		15
490901-NL	600HME, Breathing		2.0L (NL)		15
490902-NL	Breathing (2)		1.0L (NL)		20
490903-NL	Breathing (2)		2.0L (NL)		20
491000-NL	HEPA	GSL without filter	2.0L (NL)	5248	15

2x Pediatric 40"-100" (102-254 cm) Expandable Length Dual-Limb Breathing Circuit Elbow: GSP Wye: Straight				
Item Number	Filter(s)		Breathing Bag	Units per Case
490804			1.0L (L)	20
490804-NL			1.0L (NL)	20
492804	BV		1.0L (L)	20
492804-NL	BV		1.0L (NL)	20

2x Pediatric 20"-50" (51-127 cm) Expandable Length Dual-Limb Breathing Circuit Elbow: GSP Wye: Straight				
Item Number	Filter(s)		Breathing Bag	Units per Case
480804			1.0L (L)	40
482804	BV		1.0L (L)	40
482804-NL	BV		1.0L (NL)	40
482802	BV		2.0L (L)	40
482802-NL	BV		2.0L (NL)	40



Fixed Length Dual-Limb Circuits

ORDERING INFORMATION

2x Adult 40" (102 cm) Fixed Length Dual-Limb Breathing Circuit					
Elbow: GSP Wye: Straight					
Item Number	Filter(s)	Breathing Bag	Breathing Mask	Units per Case	
225-3709-804	BV	3.0L (L)		20	
370936	(2) BV	3.0L (L)		20	
370936-NL	(2) BV	3.0L (NL)		20	
670904	BV	3.0L (NL)		20	
370924	BV	3.0L (L)	5045	15	
Elbow: GSP WYE: Swivel					
370374-SG	BV	3.0L (L) SureGrip™ Bag		20	
370374-NL	BV	3.0L (NL)		20	
Elbow: Straight, GSP WYE: Swivel					
225-3703-810	BV	3.0L (L)		20	
Elbow: Standard WYE: Straight					
225-3709-800	BV	3.0L (L)		20	
Elbow: - WYE: Swivel with port					
373600	BV	3.0L (L)		20	
Elbow: - WYE: Swivel					
225-3703-800	BV	3.0L (L)		20	

ORDERING INFORMATION

2x Adult 40" (102 cm) Fixed Length Dual-Limb Breathing Circuit					
Elbow: GSP Wye: Straight					
Item Number	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
370941	BV	GSL without filter	3.0L (L)	5045	15
370941-NL	BV	GSL without filter	3.0L (NL)	5045	15
370956D	BV	GSL without filter	3.0L (L)		20
370956D-NL	BV	GSL without filter	3.0L (NL)		20
370928	BV	GSL with filter	3.0L (L)		20
370928-NL	BV	GSL with filter	3.0L (NL)		20
371056		GSL without filter	3.0L (L)		20
371056-NL		GSL without filter	3.0L (NL)		20
225-3710-804			3.0L (L)		20
671004			3.0L (NL)		20
225-3710-807			2.0L (L)		20

ORDERING INFORMATION

2x Adult 40" (102 cm) Fixed Length Dual-Limb Breathing Circuit Elbow: Standard Wye: Straight

225-3710-802	2.0L (L)	20
225-3710-800	3.0L (L)	20

2x Adult 40" (102 cm) Fixed Length Dual-Limb Breathing Circuit Elbow: Standard Wye: Swivel

374000	3.0L (L)	20
225-3704-800	3.0L (L)	20

2x Adult 60" (152 cm) Fixed Length Dual-Limb Breathing Circuit Elbow: GSP Wye: Straight

Item Number	Temp Port	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
370913-NL	YES	Breathing	GSL without filter	3.0L (NL)	5045	20
370916-NL		Breathing	GSL without filter	2.0L (NL)		20
370923-NL		Breathing (2)		3.0L (NL)		20
371804		BV		3.0L (L)		15
371804-NL		BV		3.0L (NL)		15
225-3718-806	YES	BV		3.0L (L)		15
225-3718-806-NL	YES	BV		3.0L (NL)		15
225-3719-804				3.0L (L)		15
225-3719-804-NL				3.0L (NL)		15

Elbow: Standard Wye: Straight

225-3718-800		Breathing		3.0L (L)		15
225-3719-800				3.0L (L)		15

Elbow: Standard, Straight GS Connector Wye: Straight

371821	Yes	Breathing	GSL without filter	3.0L (L)		15
--------	-----	-----------	--------------------	----------	--	----



ORDERING INFORMATION

2x Adult 72" (183 cm) Fixed Length Dual-Limb Breathing Circuit Elbow: GSP Wye: Straight					
Item Number	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
370902-NL	BV (2)		3.0L (NL)		20
370905-NL	BV	GSL without filter	3.0L (NL)	5245	20
370908-NL	BV Temp Port	GSL without filter	3.0L (NL)	5045	20
370909-NL	BV HME	GSL without filter	3.0L (L)	5045	20
225-3724-800			3.0L (L)		15
Elbow: GSP Wye: Swivel					
370300-NL	BV		3.0L (NL)		20

2x Adult 144" (366 cm) Fixed Length Dual-Limb Breathing Circuit Elbow: Standard Wye: Straight					
225-3723-800	BV		3.0L (L)		15

2x Pediatric 36" (91 cm) Fixed Length Dual-Limb Breathing Circuit Elbow: Standard Wye: Straight					
Item Number	Filter(s)		Breathing Bag		Units per Case
225-3720-800	BV		1.0L (L)		20
225-3708-802			2.0L (L)		20
225-3720-802	BV		2.0L (L)		20
Elbow: GSP Wye: Straight					
225-3708-804			1.0L (L)		20
670804			1.0L (NL)		20
672004	BV		1.0L (NL)		20
225-3720-804	BV		1.0L (L)		20
372036	(2) BV		1.0L (L)		20
370807			2.0L (L)		20
372007	BV		2.0L (L)		20

2x Pediatric 54" (137 cm) Fixed Length Dual-Limb Breathing Circuit Elbow: Standard Wye: Straight					
Item Number	Filter(s)		Breathing Bag		Units per Case
372900			1.0L (L)		15
Elbow: GSP Wye: Straight					
372904			1.0L (L)		15
375204	BV		1.0L (L)		15



Fixed Length Single-Limb Circuits

ORDERING INFORMATION

2x Adult 40" (101 cm) Fixed Length Single-Limb Breathing Circuit Elbow: GSP					
Item Number	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
384000	Insp.		3.0L (L)		20
384000-NL	Insp.		3.0L (NL)		20
384001-NL	Insp.		2.0L (NL)		20
384004	Insp.	Yes	3.0L (L)	5445	20
384004-NL	Insp.	Yes	3.0L (NL)	5445	20
384006	Insp.	Yes	3.0L (L)	5465	20
384006-NL	Insp.	Yes	3.0L (NL)	5465	20

2x Adult 60" (152 cm) Fixed Length Single-Limb Breathing Circuit Elbow: GSP					
Item Number	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
386000	Insp.		3.0L (L)		15
386000-NL	Insp.		3.0L (NL)		15
386001	Insp.		3.0L (L) SureGrip™ Bag		15
386003	Insp.		3.0L (L) SureGrip™ Bag	5445	15

2x Adult 72" (183 cm) Fixed Length Single-Limb Breathing Circuit Elbow: GSP					
Item Number	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
387901					15
387000-NL	Insp.		3.0L (NL)		15

2x Pediatric 40" (101 cm) Fixed Length Single-Limb Breathing Circuit Elbow: GSP					
Item Number	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
364001-NL	Insp.		1.0L (NL)		20
364002-NL	Insp.		2.0L (NL)		20

2x Pediatric 60" (152 cm) Fixed Length Single-Limb Breathing Circuit Elbow: GSP					
Item Number	Filter(s)	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
366001-NL	Insp.		1.0L (NL)		15
366002-NL	Insp.		2.0L (NL)		15

Manual Ventilator Transport Kits					
Item Number	Tubing	Gas Sampling Line	Breathing Bag	Breathing Mask	Units per Case
385000	Manual Ventilator Transport Kit with Dial Bleed Valve, Manometer Port with Tethered Cap without Breathing Bag				25
385100	Manual Ventilator Transport Kit with Dial Bleed Valve, Manometer Port with Tethered Cap with 3.0L Breathing Bag				20

Custom Circuits

CREATING YOUR OWN

Create your own custom order to package breathing circuit accessories and components for your specific procedures and protocols.

Choose from a wide variety of components to help streamline the ordering process (circuit tubing, wye connectors, patient elbows, breathing filters, HMEs, gas sampling lines, breathing masks, breathing bags, etc.)

Ordering custom kits, in some cases, can be more cost-effective than purchasing the components separately.

- Time savings. All the items from set-up to procedural application are in one package.
- Simple ordering. A single order for each circuit contains all the components you need, improving inventory management.
- Reduced wastage. The exact number of components is provided according to patient need and department protocol.
- Reduced stock. Custom circuits can be held by your local distributor according to forecast demand.

Components	Options
Breathing Circuits	Z-wave tubing, Square wave tubing, Expandable tubing, Adult and Pediatric single limb tubing
Breathing Masks	Premium™, Premium Plus™, Premium Soft Plus™ disposable anesthesia face masks
Elbow Connectors	90° Elbow, 90° Elbow with gas sampling port and Straight connectors
HMEs and Filters	Thermovent® HEPA Filters, HMEs
Wye Connector	Standard, Surefit, Straight, Parallel and Swivel connectors
Gas Sampling	Internal, External, Filters, Adapters
Caps	Dust, Female and Male Luer lock caps
Breathing Bags	Latex, Non-Latex, SureGrip™ Breathing Bags
Ventilator Tubes	With / Without Filters



VALUE PACK Custom Kits

CREATING YOUR OWN

- Clinicians can operate faster and more efficiently with a comprehensive kit of components.
- Operating Rooms can be set up faster while lowering turnover times.
- Value Pack means fewer items to stock, helping to make inventory more manageable, reduce re-ordering, and improve overall profitability for your healthcare facility.

Here is a list of optional components for a value Pack:

Components	Options
Breathing Circuits	Standard or Custom
Breathing Masks	Premium™, Premium Plus™, Premium Soft Plus™ disposable anesthesia face masks
Breathing Bags	Latex, Non-Latex, SureGrip™ Breathing Bags
HMEs and Filters	Thermovent® HEPA Filters, HMEs
Ventilator Tubes	With / Without Filters
ECG Pads	3 Pads / Pack, 5 Pads / Pack, Back Pad
Suction Devices	Yankauer Bulb Tip, Suction Catheter with Valve, Suction Catheter without Valve, Suction Tubing
Endotracheal Devices	Sizes 7.0-8.0 mm 10 mL Syringes Intubation Stylets
Laryngoscope Blades	Disposables in Macintosh or Miller styles
Oral Airways	Nasopharyngeal Airways I.D. 6.0-9.0 mm Berman Sizes 6.0-10.0 cm Guedel Sizes 6.0-10.0 cm
Esophageal Stethoscopes	400 Series, 12-24 french 700 Series, 12-24 french 700 Series Sub, 12-24 french
Temperature Circuits	Rectal Probe, Skin Probe, ACCU-BAR™ Temperature Strip
Blood Pressure Cuffs	Single or Double Tube
Syringes	10 mL or 20 mL
Supplies	Transpore Tape, Alcohol Prep Pads, Tongue Blades, Gloves, 3" x 3" Gauze Pads



Face Masks

ORDERING INFORMATION

Portex® Premium Soft Plus™ Face Masks with Inflation Valve			
Item Number	Face Mask Size	Hook Ring	Units per Case
005471	1		30
005452	2		30
005448	3	Clear	30
005455	4	Yellow	30
005445	5	Light Blue	30
005465	6	Dark Blue	30



Portex® Disposable Face Masks without Inflation Valve			
Item Number	Face Mask Size	Hook Ring	Units per Case
15071	1	White	50
15052	2	Pink	50
15070	3	Yellow	50
15048	4	Green	50
15055	5	Red	50
15045	6	Dark Blue	50
15266	7	Fuchsia	50



Portex® Disposable Face Masks with Inflation Valve			
Item Number	Face Mask Size	Hook Ring	Units per Case
15271	1	White	50
15252	2	Pink	50
15270	3	Yellow	50
15248	4	Green	50
15255	5	Red	50
15245	6	Dark Blue	50
15265	7	Fuchsia	50



Breathing Bags

ORDERING INFORMATION

Breathing Bags, Latex					
Item Number	Description	Size	Bushing	Units per Case	
321525	Anesthesia Breathing Bag, Natural Rubber Latex	0.5L	22.0 mm	25	
321625	Anesthesia Breathing Bag, Natural Rubber Latex	1.0L	22.0 mm	25	
321725	Anesthesia Breathing Bag, Natural Rubber Latex	2.0L	22.0 mm	25	
321825	Anesthesia Breathing Bag, Natural Rubber Latex	3.0L	22.0 mm	25	



Breathing Bags, Non-Latex					
Item Number	Description	Size	Bushing	Units per Case	
670012	Anesthesia Breathing Bag, Neoprene, Non-Latex	0.5L	22.0 mm	25	
670001	Anesthesia Breathing Bag, Neoprene, Non-Latex	1.0L	22.0 mm	25	
670002	Anesthesia Breathing Bag, Neoprene, Non-Latex	2.0L	22.0 mm	25	
670003	Anesthesia Breathing Bag, Neoprene, Non-Latex	3.0L	22.0 mm	25	



Bacterial / Viral Filter

ORDERING INFORMATION

Bacterial / Viral Filter										
Item Number	Description	Tidal VOL	Resistance (mg H ₂ O at L/min)	Moisture Output (mg H ₂ O / L@Vt)	Deadspace (mL)	Weight (gm)	Media	Filter	Filter Efficiency (%)	Units per Case
002862	Breathing Filter	n/a	1.0@30 2.0@60 3.5@90	NA	32	20	Polypropylene	Electrostatic	>99.9	40
100/585/000N	Thermovent® HEPA	150-1200	1.5@30 3.3@60 5.0@90	18.6@500 15.6@1000	45	29	Paper	Hydrophobic	>99.99999	50
002873	HEPA Filter	200-1500	1.5@30	25 @500	55	40	Paper	Hydrophobic	99.9999	50



Heat and Moisture Exchangers (HMEs)

ORDERING INFORMATION

Heat and moisture Exchangers (HMEs)								
Item Number	Description	Tidal VOL (mL)	Resistance (mg H ₂ O at L/min)	Moisture Output (mg H ₂ O / L@Vt)	Deadspace (mL)	Weight (gm)	Media	Units per Case
002815	HME, Neonatal, Straight	>10	0.9@5	32@25	3	3	Hygroscopic Foam	50
002814P	HME with Gas Port, Pediatric, Angled	>45	0.2@10	31@50	13	9	Hygroscopic Foam	50
580011	Thermovent® 600 HME, Pediatric, Straight	100-600	1.0@30 2.9@60 5.7@90	27@500	11	10	Hygroscopic Foam	50
002811P	HME with Gas Port, Rectangular, Adult, Straight	90-1500	0.98@30	32@500	76	38	Hygroscopic Foam	50
100/582/000	Thermovent® 1200 HME, Adult, Straight	600-1200	0.4@30 1.2@60 2.4@90	24@1000	32	18.6	Hygroscopic Foam	20
002813	HME with Gas Port, Adult, Small, Straight	90-1500	1.0@15	31@250	26	18	Hygroscopic Foam	50
002817P	HME with Gas Port, Adult, Large, Straight	150-1500	1.8@60	31@500	53	27	Hygroscopic Foam	50
002818P	HME with Gas Port, Adult, Small, Angled	90-1500	0.15@15	31@250	31	19	Hygroscopic Foam	50
002812	HME with Gas Port, Adult, Med., Angled	120-1500	0.22@15	>31@500	42	21	Hygroscopic Foam	50
002816	HME with Gas Port, Adult, Large, Angled	200-1500	0.3@30	31@500	66	30	Hygroscopic Foam	50
002841	HME with Flex Tube, Adult, Small, Straight	400-1000	1.0@30 3.0@60 5.0@90	30@500 27@1000	70	27	Hygroscopic Foam	20
002841P	HME with Gas Port, and Flex Tube, Adult, Small, Straight	400-1000	1.0@30 3.0@60 5.0@90	30@500 27@1000	70	27	Hygroscopic Foam	20



Filtered Heat and Moisture Exchangers (HMEFs)

ORDERING INFORMATION

filtered Heat and moisture Exchangers (HMEFs)									
Item Number	Description	Tidal VOL (mL)	Resistance (mg H ₂ O at L/min)	Moisture Output (mg H ₂ O / L@Vt)	Deadspace (mL)	Weight (gm)	Media	Filter Efficiency (%)	Units per Case
002822	Filtered HME with Gas Port, Rectangular, Adult, Straight	90-1500	0.98@30	32@500	76	38	Hygroscopic Foam	Electrostatic	50
002874	HEPA Filtered HME with Gas Port, Adult, Straight	150-1500	1.48@30	31@500	55	40	Hygroscopic Foam, Paper	Hydrophobic	50
002823P	Filtered HME with Gas Port, Small Adult, Straight	120-1500	1.3@15	>31@500	33	21	Hygroscopic Foam	Electrostatic	50
002825	Filtered HME with Gas Port, Small Adult, Angled	90-1500	1.1@15	31@250	32	18	Hygroscopic Foam	Electrostatic	50
002821	Filtered HME with Gas Port, Large Adult, Angled	200-1500	1.0@30	31@500	66	31	Hygroscopic Foam	Electrostatic	50
002865	Filtered HME, Adult, Straight	400-1000	1.6@30 3.8@60 6.5@90	27@500 24@1000	32	20	Hygroscopic Foam	Electrostatic	40
002866	Filtered HME with Flex Tube, Adult, Straight	400-1000	1.6@30 3.8@60 6.5@90	27@500 24@1000	67	34	Hygroscopic Foam	Electrostatic	20
002850P	Filtered HME with Gas Sampling Port, Adult, Straight	400-1000	1.5@30 3.0@60 5.0@90	28@500 26@1000	68	38	Hygroscopic Foam	Electrostatic	20
002851P	Filtered HME with Gas Sampling Port and Flex Tube Adult, Straight	400-1000	1.5@30 3.0@60 5.0@90	28@500 26@1000	103	52	Hygroscopic Foam	Electrostatic	20

Speciality Filters										
Item Number	Description	Tidal VOL (mL)	Resistance (mg H ₂ O at L/min)	Moisture Output (mg H ₂ O / L@Vt)	Deadspace (mL)	Weight (gm)	Media	Filter	Filter Efficiency (%)	Units per Case
002291	0.5µ Hydrophobic Disc Filter (Clear)	n/a	n/a	n/a	0.75	3.8	Gore-Tex® Membrane	Direct Interception	>99.9	200



Accessories

ORDERING INFORMATION

Gas Sampling Lines				
Item Number	Description	Luer Connection	Filter	Units per Case
225-3422-800	Gas Sampling Line, 10'x.050" ID	Male / Male		50
225-3450-800	Gas Sampling Line, 10'x.050" ID	Male / Female		50
225-3421-800	Gas Sampling Line, 10'x.050" ID	Male / Male	0.5 micron	50
225-3412-800	Gas Sampling Line, 10'x.050" ID	Male / Female	0.5 micron	50
225-3415-800	Gas Sampling Line, 10'x.050" ID	Male / Female	5.0 micron	50
225-3416-800	Gas Sampling Line, 15'x.050" ID	Male / Female	5.0 micron	50
225-3492-300	Gas Sampling Line, 10'x.047" ID	Male / Male		100
225-3445-300	Gas Sampling Line, 10'x.047" ID	Male / Female		100



Patient Connectors (Elbows)		
Item Number	Description	Units per Case
225-3507-800	(A) Patient Connector, Standard Elbow, Clear	50
225-3523-804	(B) Patient Connector, Gas Sampling Elbow with Capped Port, Clear	50
225-3524-804	(C) Patient Connector, Gas Sampling Straight with Capped Port, Clear	50



Ventilator Tubing		
Item Number	Description	Units per Case
225-3714-800	Ventilator Tube, Nonconductive, 22 mm Bushings, without Filter	30
225-3717-800	Ventilator Tube, Nonconductive, 22 mm Bushings, without Filter	30
225-3715-800	Ventilator Tube, Nonconductive, 22 mm Bushings, without Filter	20
225-3716-800	Ventilator Tube, Nonconductive, 22 mm Bushings, without Filter	20



CO ₂ Absorbents		
Item Number	Description	Units per Case
008860	Sodasorb® CO ₂ Absorbent, Canister Pak (Refill Bag) (1.64 L)	12
008870	Sodasorb® CO ₂ Absorbent, Pre Pak (Canister Discs) (1.3 L)	12
008880	Sodasorb® CO ₂ Absorbent, Pail (5 Gallon) (19 L)	1
008865	Sodasorb® LF CO ₂ Absorbent, Canister Pak (Refill Bag) (1.64 L)	12
008875	Sodasorb® LF CO ₂ Absorbent, Pre Pak (Canister Discs) (1.3 L)	12
008855	Sodasorb® LF CO ₂ Absorbent, Pail (1 Gallon) (4.3L or 4.2 kg)	3



Index

ITEM NUMBERS

100/582/000	16	002813	16	15270	13
100/585/ 000N	15	002814P	16	15271	13
225-3412-800	18	002815	16	321525	14
225-3415-800	18	002816	16	321625	14
225-3416-800	18	002817P	16	321725	14
225-3421-800	18	002818P	16	321825	14
225-3422-800	18	002821	17	364001-NL	10
225-3445-300	18	002822	17	364002-NL	10
225-3450-800	18	002823P	17	366001-NL	10
225-3492-300	18	002825	17	366002-NL	10
225-3507-500	11	002841	16	370300-NL	9
225-3507-800	18	002841P	16	370374-NL	7
225-3523-804	18	002850P	17	370374-SG	7
225-3524-804	18	002851P	17	370807	9
225-3703-800	7	002862	15	370900-NL	4
225-3703-810	7	002865	17	370902-NL	9
225-3704-800	8	002866	17	370905-NL	9
225-3708-802	9	002873	15	370908-NL	9
225-3708-804	9	002874	17	370909-NL	9
225-3709-800	7	005445	13	370913-NL	8
225-3709-804	7	005448	13	370916-NL	8
225-3710-800	8	005452	13	370923-NL	8
225-3710-802	8	005455	13	370924	7
225-3710-804	7	005465	13	370928	7
225-3710-807	7	005471	13	370928-NL	7
225-3714-800	18	008855	18	370936	7
225-3715-800	18	008860	18	370936-NL	7
225-3716-800	18	008865	18	370941	7
225-3717-800	18	008870	18	370941-NL	7
225-3718-800	8	008875	18	370956D	7
225-3718-806	8	008880	18	370956D-NL	7
225-3718-806-NL	8	15045	13	371001-NL	4
225-3719-800	8	15048	13	371002	4
225-3719-804	8	15052	13	371056	7
225-3719-804-NL	8	15055	13	371056-NL	7
225-3720-800	9	15070	13	371804	8
225-3720-802	9	15071	13	371804-NL	8
225-3720-804	9	15245	13	371821	8
225-3723-800	9	15248	13	372007	9
225-3724-800	9	15252	13	372036	9
002291	17	15255	13	372900	9
002811P	16	15265	13	372904	9
002812	16	15266	13	373600	7

ITEM NUMBERS - CONT'D

374000	8	453600-NL	5
375204	9	454000	5
384000	10	454204	5
384000-NL	10	455300	5
384001-NL	10	455417D	5
384004	10	455477D	5
384004-NL	10	480804	6
384006	10	482802	6
384006-NL	10	482802-NL	6
385000	10	482804	6
385100	10	482804-NL	6
386000	10	487000	4
386000-NL	10	490804	6
386001	10	490804-NL	6
386003	10	490900-NL	6
387000-NL	10	490901-NL	6
387901	10	490902-NL	6
450904	5	490903-NL	6
450904-SG	5	491000-NL	6
450906	5	492804	6
450906-NL	5	492804-NL	6
450907	5	580011	16
450907-NL	5	650904	5
450950-NL	4	670001	14
450951-NL	4	670002	14
450952-NL	4	670003	14
450953-NL	4	670012	14
450954-NL	4	670804	9
450955-NL	4	670904	7
450956-NL	4	671004	7
450960-NL	5	672004	9
450961-NL	5	4509901D	5
450962-NL	5		
450963-NL	5		
450964-NL	5		
450965-NL	5		
451004	5		
451004-NL	5		
451040-NL	4		
451041-NL	4		
451050-NL	5		
451051-NL	5		
453600	5		

Notes:

For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

MHYTCA-1204

Smiths Medical International
1500 Eureka Park, Lower Pemberton
Ashford, Kent, TN25 4BF
Tel: +44 (0)845 850 0445
www.smiths-medical.com

Rx
ONLY
smiths medical

Find your local contact information at: www.smiths-medical.com/customer-support/contact-us

Products referenced are CE-marked. Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Portex®, and the Smiths Medical design marks are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. ©2017 Smiths Medical. All rights reserved. SM196462GB-022017